

**Fitting Instructions:**

1. Unfasten buckles on chest support and slip neck strap over dog's head. And position bellyband between dog's legs.
2. Loop girth straps over body and fasten buckles. Adjust neck and girth straps to fit snugly.
3. Unbuckle support band and position underneath dog's hindquarters in front of rear legs and fasten buckles. End of carrying strap should be centered between buckles. Adjust to fit snugly.
4. Clip carrying strap to D-ring on chest support.
5. Lift up on carrying strap to support dog's hindquarters as it stands and walks.

**For Use as Car Restraint:**

- After fitting chest support, slide car seat belt through the nylon loop on the back strap and fasten seat belt clip.

**For Use as a Walking Harness:**

- After fitting chest support, attach lead or carrying strap to the D-ring on back.

You do not need to attach the hindquarter support band when using as a car harness.

The hindquarter support band is not usually needed when harness is being used for walking, however it can be used to provide additional support to pets with very limited mobility.

**Caution:** Recommended for pets only.

**Care:** Clean with a damp sponge and line dry.